



CRYOTHERAPY BRACE

With compression

INDICATIONS


The brace is designed for the stifle, but turned upside down, it can be used for other applications such as the elbow. Here is a non-exhaustive list of applications: elbow and stifle arthroscopy, arthrotomy, CCL surgery, fractures of the proximal tibia, fractures of the distal femur, distal femoral osteotomies to correct angular limb deformities, elbow fractures.

BENEFITS


The cold compression brace reduces inflammation, oedema and pain : it is enhanced by compression.

“COLD COMPRESSION THERAPY PROVIDES
ANALGESIA AND A FASTER RECOVERY.”

Mechanism of action



“A recent meta-analysis showed that **Cold Compression Therapy** is superior to cryotherapy alone in the acute rehabilitation stages in people undergoing knee surgery. These improvements are thought to be due to **more intensive contact of the cooling aggregate with the skin and the direct effects of compression.**”*



In dogs, when used post-operatively, the brace provides ...	24 hours after surgery	10 days after surgery	42 days after surgery
Significantly greater range of motion	✓		✓
Significantly lower Glasgow Pain Scale scores	✓		✓
Significantly improved degree of lameness		✓	✓

“The importance of rehabilitation and multimodal and pre-emptive analgesia after orthopaedic surgery has been recognised in veterinary medicine, and as such, treatment modalities without adverse effects are desirable.

Cold Compression Therapy represents such a modality, as confirmed by the present and previous studies.”*

*N. von Freeden et al. : Comparison of two cold compression therapy protocols after tibial plateau leveling osteotomy in dogs, Tierarztl Prax Ausg K Kleintiere Heimtiere 2017; 45(04): 226-233.

Comfortable for the dog and easy to apply.

Apply for 10-20 minutes, 2-3 times a day by inflating the brace with 5 to 10 pumps. For at least the first 5 days after surgery or asap after the injury.

Warning

Protect skin with a tubigrip, a gauze or bandage. Do not use with open wounds, ischemic area.

CRYOTHERAPY BRACE



Products Information

Size	Orientation	Ref.	Characteristics
S	RIGHT	AC-GDS01	15-25Kg High : 26 cm Circumference tight : 31 cm
	LEFT	AC-GGS01	
M	RIGHT	AC-GDM01	25-35Kg High : 28 cm Circumference tight : 35 cm
	LEFT	AC-GGM01	
L	RIGHT	AC-GDL01	35-50Kg High : 30 cm Circumference tight : 38 cm
	LEFT	AC-GGL01	



fourlimb.com.au



ortho@fourlimb.com.au



+61 8 9590 885

It's in your hands...

Fourlimb *surgical solutions*

