Deutsch  
English  
Español  
Français  
Italiano

marClamp® bipolare Klemmen  
marClamp® bipolar forceps  
marClamp® pinzas bipolares  
marClamp® pinces bipolaires  
marClamp® pinze bipolari

Gebrauchsanweisung  
Instructions for Use  
Instrucciones de uso  
Mode d'emploi  
Istruzioni per l'uso

**CE 0297**

**REF** 90-237-58-11  
Revision 03  
Date of Release: 2023-11

|                |    |
|----------------|----|
| Deutsch .....  | 9  |
| English .....  | 30 |
| Español .....  | 51 |
| Français ..... | 72 |
| Italiano ..... | 93 |

**Symbolerklärung**  
**Symbol Explanation**  
**Explicación de los símbolos**  
**Explications des symboles**  
**Spiegazione dei simboli**

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## Table of Contents

|       |                                                                                 |    |
|-------|---------------------------------------------------------------------------------|----|
| 1     | General .....                                                                   | 30 |
| 1.1   | Manufacturer .....                                                              | 30 |
| 1.2   | Hotline.....                                                                    | 30 |
| 1.3   | Adverse incident reporting requirement .....                                    | 30 |
| 1.4   | About this document .....                                                       | 31 |
| 1.5   | Abbreviations and Terms .....                                                   | 31 |
| 1.6   | Validity of this document .....                                                 | 31 |
| 2     | Scope of delivery .....                                                         | 31 |
| 2.1   | Inspection of the delivery for completeness and correctness.....                | 31 |
| 3     | Intended use.....                                                               | 32 |
| 3.1   | Intended Purpose.....                                                           | 32 |
| 3.2   | Indications.....                                                                | 32 |
| 3.3   | Contraindications .....                                                         | 32 |
| 3.4   | Clinical benefit.....                                                           | 32 |
| 3.5   | Possible adverse effects.....                                                   | 32 |
| 3.6   | Residual risks.....                                                             | 32 |
| 3.7   | Patient target group.....                                                       | 32 |
| 3.8   | Users.....                                                                      | 32 |
| 3.9   | Environmental conditions during use .....                                       | 33 |
| 3.10  | Restrictions on use .....                                                       | 33 |
| 3.11  | Warnings .....                                                                  | 33 |
| 4     | Operation/Use/Application .....                                                 | 36 |
| 4.1   | Storage and transport .....                                                     | 36 |
| 4.2   | Product Description .....                                                       | 36 |
| 4.3   | Combination products and accessories .....                                      | 37 |
| 4.4   | Before First Use and Before Each Further Use.....                               | 38 |
| 4.5   | Application / Intraoperative Procedure .....                                    | 38 |
| 4.6   | Connection to an electrosurgical unit.....                                      | 38 |
| 4.7   | Sealing.....                                                                    | 38 |
| 4.8   | After Completion of the Application.....                                        | 38 |
| 5     | Cleaning, disinfection, and sterilization .....                                 | 39 |
| 5.1   | Limitations and restrictions on cleaning, disinfection, and sterilization ..... | 39 |
| 5.2   | Pre-Treatment at the Application Location prior to Cleaning .....               | 40 |
| 5.3   | Preparation before Cleaning.....                                                | 40 |
| 5.4   | Cleaning and Disinfection .....                                                 | 41 |
| 5.4.1 | Manual Cleaning and Disinfection .....                                          | 41 |
| 5.4.2 | Machine Cleaning and Disinfection .....                                         | 41 |
| 5.4.3 | Validated Program Sequence .....                                                | 41 |

|       |                                                           |    |
|-------|-----------------------------------------------------------|----|
| 5.4.4 | Material Durability.....                                  | 42 |
| 5.4.5 | Drying.....                                               | 42 |
| 5.5   | Inspection, Functional Check, Maintenance.....            | 43 |
| 5.5.1 | Inspection and Functional Check.....                      | 43 |
| 5.6   | Packaging .....                                           | 44 |
| 5.7   | Sterilization .....                                       | 44 |
| 5.8   | Storage and transport .....                               | 45 |
| 5.9   | Environmental conditions for operation.....               | 45 |
| 6     | Maintenance .....                                         | 46 |
| 6.1   | General Notes .....                                       | 46 |
| 7     | Environmental information/disposal .....                  | 46 |
| 7.1   | Packaging .....                                           | 46 |
| 7.2   | Disposal .....                                            | 46 |
| 7.3   | National regulations.....                                 | 46 |
| 8     | Processing information according to DIN EN ISO 17664..... | 47 |

## 1 General

### 1.1 Manufacturer

Thank you for choosing a KLS Martin product.

This product has the CE mark, which means that it is compliant with the essential safety and performance requirements for medical devices pursuant to applicable European regulations.

We are the manufacturer of this product:

**KLS Martin SE & Co. KG**

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Tel. +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

### 1.2 Hotline

If you have questions about how to use this device or product or questions about clinical applications, please contact the product management team:

Tel: +49 7461 706-0

E-Mail: info@klsmartin.com

If you have technical questions or questions about maintenance contracts and training, please contact the Martin Service Center:

Tel: +49 7461 706-343

Fax: +49 7461 706-484

Email: service@klsmartin.com

## **NOTICE**

All packaging, and in some cases also the product itself, is marked with a batch number (LOT) and a catalogue number (REF). When making a product complaint, please include the LOT and REF.

### 1.3 Adverse incident reporting requirement

All serious incidents occurring in connection with the product shall be reported immediately to the manufacturer and to the competent authority.

## 1.4 About this document

**Possible danger to life of or serious injury to the patient, user or a third party if these instructions are not followed**

Read and observe the Instructions for Use completely. In particular, take all precautionary measures and heed all warnings.

This text refers to male, female and non-binary individuals equally, but not all pronouns have been used, solely for reasons of readability.

The electronic version of these instructions for use can be downloaded from [www.klsmartin.com](http://www.klsmartin.com).

## 1.5 Abbreviations and Terms

| Name    | Description                                                 |
|---------|-------------------------------------------------------------|
| B1, B2  | Bipolar outputs of the electrosurgical unit                 |
| CJD     | Creutzfeldt-Jakob disease                                   |
| G1 – G5 | Graduation levels for power setting                         |
| HF      | high-frequency                                              |
| HW      | Hardware release                                            |
| LTP     | Low-temperature plasma sterilization                        |
| W/D     | Washer-disinfector                                          |
| PU      | Packaging unit                                              |
| CSSD    | Central Sterile Services Department (central sterilization) |

## 1.6 Validity of this document

This document applies to marClamp bipolar forceps.

## 2 Scope of delivery

marClamp bipolar forceps

Catalogue numbers: 80-975-16-04, 80-975-18-04, 80-975-23-04, 80-975-28-04, 80-976-23-04, 80-977-23-04, 80-978-23-04

### 2.1 Inspection of the delivery for completeness and correctness

Immediately after receiving the delivery, check that it is complete and in good condition.

Report any damage that occurs during transit promptly.

After delivery of the product, check the original packaging and the packaging seal for intactness.

### **3 Intended use**

#### **3.1 Intended Purpose**

The bipolar sealing system marClamp is intended for permanent occlusion of veins, arteries and tissue bundles. It is designed for use in open surgery applications.

#### **3.2 Indications**

The indications result from the intended purpose.

#### **3.3 Contraindications**

- Note that tubal ligatures with the marClamp cannot be considered a reliable method for contraception.
- Not for use on the appendix!

#### **3.4 Clinical benefit**

Sealing of vessels and tissue bundles

#### **3.5 Possible adverse effects**

In general, all surgical procedures are associated with the risk of possible adverse effects and complications. In most cases, these are not caused by the product used, but rather due to clinical conditions.

There are no known product-specific adverse effects when used as intended. The user is responsible for informing the patient about possible adverse effects and complications associated with the surgical intervention.

#### **3.6 Residual risks**

The user is responsible for informing patients about the residual risks associated with the use of this product.

#### **3.7 Patient target group**

There are no restrictions regarding the patient target group.

#### **3.8 Users**

The products are to be used exclusively by qualified surgical staff.

Cleaning, disinfection and sterilization are carried out by trained qualified staff in the processing unit for medical products.

### 3.9 Environmental conditions during use

Use only in ORs under conditions appropriate for surgery or in medical facilities designed for this.

### 3.10 Restrictions on use

#### **WARNING**

##### **Danger of injury in case of insufficient qualification of the user!**

Safe use of the marClamp in combination with the current type SealSafe® IQ requires the user to be familiar with the technology and applications of electrosurgery.

Instruments for electrosurgery may be used only by persons who have been specially trained or instructed in their use.

#### **CAUTION**

##### **Danger of injury from improper application and abuse or misuse!**

Note that tubal ligatures with the marClamp cannot be considered a reliable method for contraception.

Not for use on the appendix.

### 3.11 Warnings

#### **WARNING**

##### **Non-compliance with the safety measures may result in serious or even fatal injury to the patient!**

#### **WARNING**

##### **Danger of explosion and fire from flammable gases!**

Never use flammable anesthetics or oxidizing gases when working with the marClamp bipolar sealing system.

- Before sealing, ensure that no endogenous gases are present in the area of application.

#### **WARNING**

##### **Risk of burns upon contact of the active electrode with metallic parts!**

If the active electrode is allowed to get in contact with metal parts, shunts or concentrated leakage paths can be created for the HF current. These can cause burns.

- Do not touch any other metallic instruments or objects during activation.

 **WARNING****Danger of electromagnetic interference in the presence of active implants!**

In patients with pacemakers or other active implants, there is danger of interference with or damage to the active implant.

- Before performing the surgery, consult a cardiologist and the manufacturer of the pacemaker or active implant.

 **WARNING****Danger of injury from hot surfaces!**

The electrodes in the jaw part and the outer surfaces of the jaw part are so hot due to the heat conduction even after the sealing procedure that they can cause unintentional burns.

- Avoid contact with the jaw part and tissue contact immediately after sealing.

 **WARNING****Danger of injury from unintentional contact with the patient's tissue!**

Unintentional activation via the foot switch can lead to burns and electric shocks if the electrosurgical device is placed on the patient and thus comes into contact with tissue.

- Do not place temporarily unused instruments on the patient in the intervals between use.
  - Place instruments separately from the patient, e.g. on the instrument table.
- Seal only if the jaw part of the instrument is located within sight! Unintentional activation or movement of the instrument may result in injury to the patient.
- Do not touch any other metallic instruments or objects during the sealing process.
- Always avoid accidental tissue contact or contact with similar low-resistance materials (e.g. liquids) of the instruments.

---

 **CAUTION**

**Danger of injury for patient and user from improper application!**

Improper application and abuse or misuse may result in injury to patient and user and/or premature wear of the sealing system.

- Clean and sterilize the marClamp before each use, see section 5, “Cleaning, disinfection, and sterilization”, page 39.
  - Use the bipolar sealing system only for its intended purpose, see section 3.1, “Intended Purpose”, page 32.
  - Route the connection cable so that it touches neither the patient nor other cables, and does not form an obstacle. Protect the cable from mechanical damage (running over, crushing, bending).
  - Use of high-frequency devices can interfere with other devices.
- 

---

 **WARNING**

**Possible danger to life for third parties in case of shipment of contaminated products!**

When returning products, send only cleaned and disinfected products in sterile packaging.

---


## 4 Operation/Use/Application


### 4.1 Storage and transport

- Store the bipolar forceps in a clean, cool and dry place.
- Protect from mechanical damage.
- Store and transport in safe containers/packing. Use certified sterilization packing (e.g. according to EN 868, ISO 11607) or sterilization containers for sterilization. We recommend the KLS Martin screen basket system.
- Handle with great care. Do not drop or throw.
- Keep away from sunlight.
- In case of return shipments send only clean and disinfected products in sterile packing.

### 4.2 Product Description

marClamp bipolar forceps are available in varying lengths and configurations. Figures can be found below.

### 4.3 Combination products and accessories

#### **NOTICE**

The bipolar sealing instrument marClamp is intended exclusively for use with the electrosurgical unit maxium® (REF 80-042-00-04, REF 80-042-02-04 and REF 80-042-04-04, type ME 402) from software version V3.392, or new developments by the manufacturer. The IQ connector of the marClamp Cut can be connected to the bipolar socket B1 or B2 of the electrosurgical units maxium®. The current type SealSafe® IQ is selected automatically. It can be used exclusively with IQ instruments (Plug & Play)!

#### **NOTICE**

For correct operation and adjustment, observe the Instructions for Use of the electrosurgical unit maXium®.

- Connect connection cable of the marClamp to the bipolar HF output B1 or B2 of the electrosurgical unit maxium®.
  - The electrosurgical unit maXium® automatically supplies the required current type SealSafe® IQ with a factory default setting G3.

The power setting is adjusted in graduation levels (G grades). It can be varied from G1 to G5. Here, with increasing G grade the switch-off impedance is increased and the sealing time extended. The foot switch for the activation of the sealing current must be selected separately.

---

#### 4.4 Before First Use and Before Each Further Use

- Observe the safety notes, see section 3.1, “Intended Purpose”, page 32.
- Perform a functional check before each use (e.g. on a swab soaked in NaCl solution).
- Inspect the instrument for damage or missing components prior to each use (the protective cap over the instrument joint may come loose if cleaned improperly).

---

 **WARNING**

##### **Risk of infection in case of non-sterile handling!**

Non-sterile handling and improper sterilization can cause severe health risks for the patient.

The instruments must be cleaned and sterilized before the first and any further use.

---

#### 4.5 Application / Intraoperative Procedure

Before each use, a visual inspection and functional check must be performed on the forceps, after which they must be cleaned and sterilized. Inspect the bipolar forceps by checking for cleanliness, intact insulation, and integrity before each use. Withdraw damaged forceps from service immediately!

The coagulation current is activated using a foot switch.

#### 4.6 Connection to an electrosurgical unit

---

**NOTICE**

Please follow the instructions for use for the relevant electrosurgical unit for proper use / configuration.

---

#### 4.7 Sealing

- Grasp the tissue with the instrument tips.
- Keep in mind that the tissue structures to be sealed must not be too thin. For the sealing process, the amount of tissue must sufficiently fill the space between the branches after the lock has latched!
- The tips of the forceps must remain within the surgeon's field of view.

#### 4.8 After Completion of the Application

- Disconnect the bipolar cable from the electrosurgical unit.

## 5 Cleaning, disinfection, and sterilization

### **WARNING**

#### **Possible danger to life of the patient in case of non-sterile handling!**

Clean, disinfect and sterilize products supplied in non-sterile condition before the first and before any further use.

- In case of patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible CJD variants, processing of the products must be carried out according to the applicable national regulations.
- The operator/processor is responsible for cleaning, disinfecting and sterilizing the products used. It is essential that national regulations, including restrictions, be observed.

The processing procedure must be validated prior to use of the medical device. The operator/processor is responsible for this.

Because of potentially inadequate labeling on the product, products that have been removed from their sterile packaging but have not been used and are not contaminated cannot be considered equivalent to products that have been obtained as non-sterile products. In this case, the operator/processor is responsible for use including labeling, cleaning, disinfection and sterilization.

The specific information for each reference number for the following section is available at:

[www.klsmartin.com/processing](http://www.klsmartin.com/processing)

### **5.1 Limitations and restrictions on cleaning, disinfection, and sterilization**

The instrument is not subject to any restrictions on the number of reprocessing cycles. With each application, the wear of the product progresses. Product lifetime is usually determined by wear and damage during use. This wear is due to technical reasons and is unavoidable. Before each use, a visual inspection and functional check must therefore be performed on the forceps.

Based on practical experience and user observations, the products can be reprocessed 50 times. However, the instrument is typically used far past this number, with the actual number of uses highly dependent on the degree of care in treatment by the user's different departments.

We would like to note here that the instrument must be handled and reprocessed as instructed in the applicable instructions for use.

**NOTICE****Risk of damaging the forceps!**

- Do not immerse the bipolar forceps into hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>).
- To avoid scratching the surfaces of the tips and the insulation, do not use abrasive cleaning tools such as metal brushes or steel wool.

**5.2 Pre-Treatment at the Application Location prior to Cleaning**

We generally recommend dry return transport in a closed transport container without the addition of liquids, cleaning agents or disinfectants.

**5.3 Preparation before Cleaning**

The operator or the product reprocessor are responsible for ensuring the proper outcome of cleaning, disinfection, and sterilization of the instruments used.

Ensure that

- only product-specific and adequately validated cleaning, disinfection, and sterilization procedures are used,
- the equipment used (disinfector, sterilizer) is serviced and inspected at regular intervals and that validated parameters are used in each cycle,
- National guidelines and regulations as well as the hygiene rules of the medical practice or hospital are observed (especially with regard to prion inactivation).

When selecting the cleaning agent system, it must be ensured that

- it is suitable for cleaning instruments made of metals and plastics,
- a suitable disinfectant with verified effectiveness (e.g. VAH/DGHM or FDA approval or CE marking) is used and that it is compatible with the cleaning agent used, and
- the chemicals used are compatible with the instruments (see section 5.4.4, “Material Durability”, page 42).

Clean bipolar forceps thoroughly under running water or with a aldehyde-free disinfectant directly after use (within a maximum of 2 hours) to prevent particles and dried secretions from sticking to the instruments. Blood and tissue residues should not be allowed to dry on the instruments.

Due to the better results in terms of efficacy and repeatability, machine cleaning should be preferred to manual cleaning.

## 5.4 Cleaning and Disinfection

### 5.4.1 Manual Cleaning and Disinfection

Recommendations for manual cleaning: Validated cleaning and disinfection procedure acc. to DIN EN ISO 15883. Deviating cleaning procedures, validated acc. to DIN ISO 15883, are possible.

Manual pre-cleaning serves as preparation for machine cleaning.

### 5.4.2 Machine Cleaning and Disinfection

Recommendations for machine cleaning: Validation acc. to DIN EN ISO 15883 and information to be provided acc. to DIN EN ISO 17664. Cleaning with other process parameters, validated according to DIN EN ISO 15883, is possible.

## **NOTICE**

### **Risk of damaging the forceps!**

When selecting the disinfectant, it must be ensured that

- it has a verified effectiveness (e.g. DGHM or FDA approval or CE marking in accordance with DIN EN ISO 15883),
- the program used is suitable for the instruments and includes a sufficient number of flushing cycles,
- only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water) is used for rinsing,
- the air used for drying is filtered, and
- the disinfectant is serviced and inspected at regular intervals.

Adhere to the criteria specified in chapter 8, "Processing information according to DIN EN ISO 17664", page 47 when selecting the cleaning agent system to be used.

Place the instruments into the disinfectant. Take care to ensure that the instruments do not touch each other.

### 5.4.3 Validated Program Sequence

#### **Pre-cleaning**

**Rinsing temperature:  $10 \pm 2$  °C**

**Insertion time: 1 min**

#### **Cleaning**

Cleaning temperature:  $70 \pm 2$  °C

Insertion time: 5 min

Cleaning agent: deconex ALKA One-x, concentration 0.5 % (5 ml/l)

**Post-rinsing**

Rinsing temperature:  $10 \pm 2$  °C

Insertion time: 1 min

**Thermal disinfection**

Disinfection temperature:  $90 \pm 2$  °C

Insertion time: 5 min

Inspect the instruments as soon as possible after removal for visible contamination and package after an additional drying cycle, whenever necessary, at a clean site.

**5.4.4 Material Durability**

Make sure that the cleaning or disinfection agent does not comprise any of the following substances:

- Organic acids, mineral acids, oxidizing acids (lowest allowable pH value 5.5)
- Strong lye (highest allowable pH value 12)
- Organic solvents (e.g. acetone, ether, alcohol, gasoline)
- Oxidants (e.g. peroxides)
- halogens (chloride, iodine, bromine)
- Aromatic, halogenated hydrocarbons

**5.4.5 Drying**

---

**NOTICE****Ensure reliable drying!**

Reliable drying is crucial for the success of subsequent sterilization.

---

## 5.5 Inspection, Functional Check, Maintenance

### 5.5.1 Inspection and Functional Check

#### **NOTICE**

Whether or not the product can be used any further shall be confirmed by successful inspection of the product. The inspection release and packing into a sterile barrier system releases the product for the next application.

- The products must be macroscopically clean, i.e. free from visible contamination, after each cleaning process.
- Stop using any instruments with stains immediately and subject them to special treatment.
- In case of defects or damage, discard the instruments immediately.
- We recommend storing the bipolar forceps in a mesh basket after the inspection to prevent transportation damage and using suitable protective devices to protect them from damage.
- Perform a visual check for intact insulation, cleanliness and potential damage of the forceps.

## 5.6 Packaging

Approved sterilization packaging (e.g. conforming to EN 868, ISO 11607) must be used for sterilization, subsequent transportation and storage.

## 5.7 Sterilization

### **WARNING**

#### **Peril to life or danger of serious physical injury due to prion infection!**

Destroy bipolar forceps immediately after potential contact with prions (risk of CJD contamination!) and do not use again under any circumstances.

### **WARNING**

#### **Risk of infection in case of non-sterile handling!**

Improper sterilization and non-sterile handling of the forceps can lead to serious health hazards for patients.

Sterilization must be carried out according to a validated steam sterilization process, for example in a sterilizer complying with EN 285:2009 and DIN EN 13060, validated in accordance with ISO 17665-1:2006.

### **NOTICE**

#### **Risk of damaging the forceps!**

- Do not sterilize bipolar forceps using hot air.
  - Do not sterilize bipolar forceps in a “STERRAD”.  
(STERRAD is a registered trademark of Johnson & Johnson Inc.).
  - Do not use “flash” autoclaving.
  - Do not use irradiation sterilization, gas sterilization, or formaldehyde sterilization.
- 
- Use certified sterilization packing (e.g. according to EN 868, ISO 11607) or sterilization containers for sterilization. We recommend the KLS Martin screen basket system.
  - Sterilization at high temperatures and prolonged sterilization times shortens the life span of the forceps.
  - Ensure sufficient drying.
  - As regards sterilizer handling, loading and parameter setting, please follow the recommendations of the manufacturer of your machine.
  - The operator must ensure that sterility is maintained after the sterilization process.

Only use the sterilization methods listed below. Other sterilization procedures are not permitted.

| Method                                       | Exposure temperature       | Exposure time             | Drying time |
|----------------------------------------------|----------------------------|---------------------------|-------------|
| Fractionated vacuum<br>(dynamic air removal) | min. 132 °C<br>max. 138 °C | min. 3 min<br>max. 20 min | min. 30 min |
| Gravity method <sup>1</sup>                  | min. 132 °C                | min. 40 min               | min. 30 min |

- Steam sterilizer according to EN 13060 and EN 285
- validated according to DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance evaluation (PQ))
- make sure not to exceed a sterilization temperature of 138 °C (plus tolerance according to DIN EN ISO 17665)

**NOTE:** The fundamental suitability of the instruments for effective steam sterilization by the fractionated vacuum procedure and the gravity procedure was demonstrated by an independent accredited test laboratory. Typical conditions in clinics and physicians' practices as well as the procedure described above were taken into account.

## 5.8 Storage and transport

Information on minimum and maximum temperature for storage and transportation.

| Name                                   | marSeal 5 plus                                 |                      |
|----------------------------------------|------------------------------------------------|----------------------|
| Environmental conditions for storage   | Ambient temperature                            | +4.4 °C up to +49 °C |
|                                        | Relative humidity<br>(no condensing moisture!) | 10 – 70%             |
| Environmental conditions for transport | Ambient temperature                            | –27 °C up to +58 °C  |
|                                        | Relative humidity<br>(no condensing moisture!) | 35 – 80%             |

- Store the instruments in a clean, cool and dry place.
- Keep away from sunlight.
- Protect instruments from mechanical damage.
- Store and transport instruments in safe containers/packing. Appropriately approved sterilization packaging (e.g. according to DIN EN 868, ISO 11607) or sterilization containers must be used for this purpose.
- Handle instruments with great care. Do not drop or throw.

## 5.9 Environmental conditions for operation

See IFU HF generator maXium.

<sup>1</sup> The less-effective gravity method may only be used if the fractionated vacuum method is not possible.

## **6 Maintenance**

### **6.1 General Notes**

The product may be repaired only by the manufacturer or a qualified person or firm expressly authorized by the manufacturer to perform such work.

Modification of the product can lead to unforeseeable risks and is therefore not permitted.

If the repair is carried out by a person or firm specially authorized by the manufacturer, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must show the date of the repair and the details of the person or firm carrying out the work and must be signed.

In all cases where a party other than the manufacturer performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product are never allowed and shall void any claims of liability against the manufacturer.

## **7 Environmental information/disposal**

### **7.1 Packaging**

The manufacturer will take back all of the packaging on request. Whenever possible, parts of the packaging will be reused.

If this is not an option, the packaging can be disposed of with paper/household waste.

### **7.2 Disposal**

The product was designed to use as few composite materials as possible. This design approach ensures that it can be recycled to a large degree. This is why we offer to take back the product and to dispose of it properly.

### **7.3 National regulations**

Follow national waste management regulations and guidelines for all disposal activities.

## 8 Processing information according to DIN EN ISO 17664

|                                                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <b>WARNING</b> | <p><b>Danger of prion infection and risk of damage to the bipolar forceps from improper use!</b></p> <ul style="list-style-type: none"> <li>• <b>Destroy forceps immediately after potential contact with prions (risk of CJD contamination) and do not use them again!</b></li> <li>• <b>Do not use the “STERRAD” plasma sterilization method!</b></li> <li>• <b>Do not sterilize with hot air!</b></li> <li>• <b>Do not use “flash” autoclaving!</b></li> <li>• <b>Do not immerse the instruments into hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>)!</b></li> <li>• <b>To avoid scratching polished working tips and insulation, do not use abrasive cleaning tools!</b></li> </ul> |
| <p><b>Reprocessing limitations</b></p>                                                           | <p>Damage such as scratched contact surfaces, loss of insulation due to cracks or splintered plastic coating and wear.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

### INSTRUCTIONS

|                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Workplace</b></p>                                     | <p><b>CSSD, instrument processing department</b></p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <p><b>Storage and transport</b></p>                         | <p>Store the forceps in a clean, cool and dry place. Protect from mechanical damage. Handle with great care. Do not drop or throw.<br/>Store and transport in safe containers/packing.</p>                                                                                                                                                                                                                                                                                                                                                                                                  |
| <p><b>Preparation for cleaning</b></p>                      | <p>Clean forceps thoroughly under running water or with a aldehyde-free disinfectant directly after use (within a maximum of 2 hours) to prevent particles and dried secretions from sticking to the instruments.</p>                                                                                                                                                                                                                                                                                                                                                                       |
| <p><b>Cleaning: Manually</b></p>                            | <p>There is no specified validated manual procedure, so the machine procedure should be given preference. If a mechanical procedure is not available, when selecting the cleaning agent system ensure that</p> <ul style="list-style-type: none"> <li>• it is suitable for cleaning instruments made of metals and plastics,</li> <li>• a suitable disinfectant with verified effectiveness (e.g. VAH/DGHM or FDA approval or CE marking) is used and that it is compatible with the cleaning agent used, and</li> <li>• The chemicals used are compatible with the instruments.</li> </ul> |
| <p><b>Cleaning: Mechanical procedure</b></p>                | <p><b>Validated Program Sequence:</b></p> <p>Pre-cleaning: Rinsing temperature: 10 ± 2 °C, Insertion time: 1 min<br/>Cleaning: Cleaning temperature: 70 ± 2 °C, Insertion time: 5 min, Cleaning agent: deconex ALKA One-x, concentration 0.5 % (5 ml/l)<br/>Final rinse: Rinsing temperature 10 ± 2 °C, Insertion time: 1 min</p>                                                                                                                                                                                                                                                           |
| <p><b>Disinfection</b></p>                                  | <p>Thermal disinfection for mechanical treatment 90 ± 2 °C, insertion time 5 minutes.<br/>Immersion in disinfectant solution. Observe manufacturer's information on concentration, exposure time and temperature. Observe material compatibility (see section 5.4.4, page 42), 90°C, min. 5 minutes.</p>                                                                                                                                                                                                                                                                                    |
| <p><b>Maintenance, checks and functional inspection</b></p> | <p>Check forceps for damage prior to each use.<br/>Withdraw damaged forceps from service immediately.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

|               |                                                                                                                                       |                            |                           |             |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------|----------------------------|---------------------------|-------------|
| Packaging     | Use certified sterilization packing (e.g. according to EN 868, ISO 11607) or sterilization containers.                                |                            |                           |             |
| Sterilization | Prior to sterilization, the forceps must be cleaned, disinfected, rinsed with distilled water and finally dried carefully.            |                            |                           |             |
|               | Method                                                                                                                                | Exposure temperature       | Exposure time             | Drying time |
|               | Fractionated vacuum (dynamic air removal)                                                                                             | min. 132 °C<br>max. 138 °C | min. 3 min<br>max. 20 min | min. 30 min |
|               | Gravity method                                                                                                                        | min. 132 °C                | min. 40 min               | min. 30 min |
|               | As regards sterilizer handling, loading and parameter setting, please follow the recommendations of the manufacturer of your machine. |                            |                           |             |
| As of 03/2020 |                                                                                                                                       |                            |                           |             |

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