


IXOBONE[®]

BIOMATERIALS

BIOMIMETIC | OSTEOSTIMULATIVE | STRENGTH | READY TO USE | EXCELLENT HANDLING



Bone Composition




The two major components of bone are:

- Calcium phosphate
- Collagen


Four types of cell are present in bone:

- Osteoblasts – synthesise and deposit new bone matrix
- Osteocytes – maintains bone tissue
- Osteoclasts - resorb bone by direct chemical and enzymatic attack
- Undifferentiated stem cells transform into osteoblasts under appropriate stimuli

Bone Repair




Stage 1: Inflammatory molecules present in the blood clot recruit cells essential for the repair process.



Stage 2: Connective tissue infiltrates into the defect forming a bridge between the bone surfaces. New blood vessels form and penetrate the defect.



Stage 3: Osteoblasts form woven bone fully bridging the gap.



Stage 4: The new bone is remodelled replacing woven bone with cortical bone.

IXOBONE biomaterials can be mixed with bone marrow aspirate or autologous bone.

IXOBONE biomaterials act as a bridge between the bone surface, encouraging rapid cell colonisation and angiogenesis

IXOBONE biomaterials act as scaffolds both promoting and supporting new bone formation.

IXOBONE biomaterials will be resorbed or remodelled overtime.

IXOBONE® Product Range



IXOBONE® HA/TCP BLOCKS/DISCS/WEDGES/GRANULES

COMPOSITION:	60% HA, 40% TCP
POROSITY:	81%
MACROPOROSITY:	200 µm – 800 µm
FULLY INTERCONNECTED:	YES
MICROPOROSITY:	1.0 µm – 10 µm
STRENGTH:	5.5 to 7.5 Mpa
STERILITY:	Gamma irradiated

IXOBONE® PASTE

HIGH SURFACE AREA:	Approx. 100 m ² /g
HYDROXYAPATITE CONTENT:	38%
HA NANOTECHNOLOGY:	30-50 nm particles
PRODUCT VOLUME:	0.5 cc to 15 cc
STERILITY:	Gamma irradiated

IXOBONE® PUTTY

HYDROXYAPATITE CONTENT:	38%
HIGH SURFACE AREA:	100 m ² /g
HA NANOTECHNOLOGY:	30-50 nm particles
MICRO-SCAFFOLDS CONTENT:	16.5%
GRANULAR MICRO-SCAFFOLDS:	0.5-1.0 mm
PRODUCT VOLUME:	0.5 cc to 15 cc
STERILITY:	Gamma irradiated


IXOBONE[®] Product Range

Technical Overview



IXOBONE[®] HA/TCP

- High Porosity
- Physico-chemical characteristics tailored for maximum benefit
- Interconnected pores
- Excellent Integral Strength
- Resorbable
- Osteoconductive



IXOBONE® HA/TCP

Indications for Use




- Similarity of the internal structure to that of cancellous bone
- Exceptional porosity versus strength ratio due to proprietary technology
- Pre-shaped in standard sizes, sterile and easy to handle
- Greater interconnected structure than in other synthetic products
- Both macro and micro porosity
- High integral strength



IXOBONE® HA/TCP – Product Benefits

VET

- Easy to use
- Blocks and wedges can be easily cut and shaped to fit the defect
- Can be soaked with blood or bone marrow aspirate
- Granules can be mixed with autologous bone
- No risk of disease transmission
- Synthetic - available in unlimited quantities
- May eliminate the need for the harvesting of autologous bone









SEM showing IXOBONE® HA/TCP
soaked with blood



 **exabone**

IXOBONE[®] HA/TCP – Instructions for Use

VET



Step 1: Remove product from packaging

Step 2 (optional): The implant may be mixed with saline, blood, bone marrow aspirate or autograft

Step 3: Ensure product is fully soaked

Step 4: Fill the defect with the bone graft ensuring contact with existing healthy bone. Close the wound



 **exabone**




IXOBONE® HA/TCP – Contraindications include:

- Must not be used to provide structural support or to gain screw fixation
- Existing acute or chronic infections
- Severe vascular, neurological, degenerative disease or uncontrolled diabetes
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease or malignant tumours
- Severely impaired renal function



IXOBONE[®] PASTE

- 30 - 50nm hydroxyapatite particles
- Ultra-high surface area (~100M²/g)
- Easy to use injectable formulation
- Readily resorbable
- Osteostimulative



IXOBONE[®] PASTE

Indications for Use



- Mouldable/injectable
- Sticky consistency maximises the bone-implant interface
- Ultra-high surface area – adsorbs the biomolecules essential for bone repair
- Osteostimulative effect
- Highly resorptive nature results in locally elevated levels of Ca^+ and PO_4^- ions encouraging bone regeneration



IXOBONE® PASTE - Product Benefits

- Ready to use
- Delivery device allows for application directly into the defect
- Can be mixed with bone marrow aspirate and/or autologous bone
- No risk of disease transmission
- Available in unlimited quantities
- May eliminate the need for the harvesting of autologous bone and donor site morbidity



IXOBONE[®] PASTE – Instructions for Use

VET



Step 1: Remove product from packaging and unscrew cap



Step 2: Securely screw the most appropriate cannula onto syringe



Step 3: Inject product directly into the defect site



Step 4: Remove any remaining product from the cannula using the clearance tool



Step 5: Ensure the defect is fully filled. Close the wound



 **exabone**


IXOBONE® PASTE – Contraindications include:

- Must not be used to provide structural support or to gain screw fixation
- Existing acute or chronic infections
- Severe vascular, neurological, degenerative disease or uncontrolled diabetes
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease or malignant tumours
- Severely impaired renal function



IXOBONE[®] PUTTY

- Easy to use injectable putty
- Positive osteostimulative effect
- Multi-phasic regenerative activity
- Stratified resorption provides longer term support during bone formation and remodeling phases



IXOBONE[®] PUTTY

Indications for Use



IXOBONE® PUTTY - Features

VET

- Synthetic mix of hydroxyapatite/tri-calcium phosphate granules in an osteostimulative carrier
- Delivery device allows for application directly into the defect
- Non-setting formula
- Sticky formula provides excellent contact with host bone
- Stratified resorption at a controlled rate



 exabone

IXOBONE® PUTTY - Product Benefits

VET

- Ready to use – putty-like consistency for injection directly into defect site
- Additive free
- Synthetic, sterile and reliable alternative to autograft or allograft
- Resorption rate providing structural support throughout the regenerative process and maintaining bone volume
- Can be mixed with bone marrow aspirate and/or autologous bone



 **exabone**

IXOBONE® PUTTY – Instructions for Use

VET



Step 1: Remove product from packaging and unscrew cap



Step 2: Securely screw the most appropriate cannula onto syringe



Step 3: Inject product directly into the defect site



Step 4: Remove any remaining product from the cannula using the clearance tool



Step 5: Ensure the defect is fully filled. Close the wound



 **exabone**

IXOBONE® PUTTY – Contraindications include:

VET

- Must not be used to provide structural support or to gain screw fixation
- Existing acute or chronic infections
- Severe vascular, neurological, degenerative disease or uncontrolled diabetes
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease or malignant tumours
- Severely impaired renal function



 **exabone**

Exabone GmbH

Route de Saint-Cergue 14

CH-1260 Nyon

Switzerland

T +41-22 362 03 94

F +41-22 362 03 93

E info@exabone.com W www.exabone.com



TRADEMARKS:

- IXOBONE® is a registered trademark of Exabone GmbH in the E.U. and Switzerland, and an unregistered trademark in other territories.
- “The Swiss Bone Architect”™ and the exabone logo are unregistered trademarks of Exabone GmbH.