

BOSS Lights

Illuminate your procedures with shadowless, cool, bright light

Neuromed is delighted to introduce two new disposable lighting systems into the veterinary market

LightMat®

LightMat™ The patented fiber-optic technology provides white, shadow less light throughout the procedure, allowing the surgeon to illuminate the surgical cavity with cool, bright light.

The Lightmat range incorporates different sizes and models, and allows each device to be used as a retractor, or to be attached to an existing retractor or flat surfaced device.

VersaLight™

VersaLight™ This four-in-one hand tool increases clinical efficiency, providing not just cool white light directly where you need it, but also suction and irrigation, regulated with controls on the handle, and moderate retraction with two available handle angles.



LightMat®


Ultra-Thin LightMat w/adhesive backing

Code	Price
OMUA2550	\$731.00

Adhesive Strip:
½" x 7" (13mm x 178mm)

Illumination Area:
½" x 2" (13mm x 51mm)

Cable Size:
4mm



BOSS
INSTRUMENTS LTD

Ultra-Thin LightMat w/built-in malleable retractor

Code	Price
OMTA2334	\$803.00

Malleable Strip:
½" x 7" (13mm x 178mm)

Illumination Area:
½" x 2" (13mm x 51mm)

Cable Size:
4mm



BOSS
INSTRUMENTS LTD

LightMat®

OrthoLight Spotlight w/adhesive backing

Code	Price
OMMM1462	\$772.00

Malleable Strip

¼" x 7" (6.5mm x 178mm)

Illumination Area

¼" x 1¼" (6.5mm x 32mm)

Cable Size

2.5mm



BOSS
INSTRUMENTS LTD

OrthoLight Spotlight w/ultra-malleable backing and adhesive strip

Code	Price
OMOM2130	\$803.00

Adhesive Strip:

½" x 7" (13mm x 178mm)

Illumination Area:

½" x ½" (13mm x 51mm)

Cable Size:

2.5mm



BOSS
INSTRUMENTS LTD

VersaLight™

VersaLight Standard

Code	Price 5pack
OMVSG211 (without tubing)	\$1047.00
OMVSG311 (with tubing)	\$1209.00



Size
13½" x ½" (342.9mm x 12.7mm)

Handle Angle
15°

Cable Size
2.5mm

BOSS
INSTRUMENTS LTD

VersaLight Standard with Frazier Tip

Code	Price 5pack
OMVSX333	\$1209.00



Size
15¼" x ½" (387.4mm x 12.7mm)

Handle Angle
15°

Cable Size
2.5mm

BOSS
INSTRUMENTS LTD

VersaLight™

VersaLight Mini

Code	Price 5pack
OMVSG411	\$1161.00



Size
11¼" x ½" (285.8mm x 12.7mm)

Handle Angle
30°

Cable Size
2.5mm

BOSS
INSTRUMENTS LTD

VersaLight Mini with Frazier Tip

Code	Price 5pack
OMVSX433	\$1161.00



Size
13" x ½" (330.2mm x 12.7mm)

Handle Angle
30°

Cable Size
2.5mm

BOSS
INSTRUMENTS LTD

Light Cables

ACMI

Size	Code	Price
2.5mm	OM006241	\$850.00
2.5mm Yellow	OM007986	\$850.00
4mm	OM005011	\$850.00



Storz

Size	Code	Price
2.5mm	OM008217	\$850.00
4mm	OM005366	\$850.00

Wolf (Stryker)

Size	Code	Price
2.5mm	0OM08218	\$850.00
4mm	0OM05757	\$850.00

Olympus

Size	Code	Price
2.5mm	OM008219	\$850.00
4mm	OM005365	\$850.00

BRICK Light Source

Watts	Code	Price
50w	OM008554	\$5750.00

50w BRICK light source to power VersaLights and tubing systems.

