

## LINEAR CUTTER - GIA

Product Code			Stapling Line			Individual Staple	
Size	Stapler	Reload	Rows	Length	Number	Closed Size	Leg Length
60mm Blue	LC60B	LCR60B	4	61mm	60	1.5mm	3.8mm
60mm Green	LC60G	LCR60G	4	61mm	60	2.0mm	4.5mm
80mm Blue	LC80B	LCR80B	4	81mm	80	1.5mm	3.8mm
80mm Green	LC80G	LCR80G	4	81mm	80	2.0mm	4.5mm