




PEDICLE SCREW 5.5mm SYSTEM

3.5mm SCREW (BLUE) DRILL SIZE 2.5mm




4.5mm SCREW (GOLD) DRILL SIZE 3.2mm





Always check scale

PEDICLE SCREW 3.0mm SYSTEM




2.0mm SCREW (GOLD) DRILL SIZE 1.5mm



2.7mm SCREW (BLUE) DRILL SIZE 2.0mm



3.5mm SCREW (GREEN) DRILL SIZE 2.5mm



Always check scale