


# Elevation TPLO Plate

Locking - Polyaxial

## Correct order of screw placement



## Weight Plate Guidelines

NB Guidelines only.

Pre-operative planning with consideration of radiographic findings, conformation, patient factors and owner's compliance is essential.



Dimensions	Patient Weight
2.0mm TPLO Plate - 28mm	3.5 - 7kg
2.4mm TPLO Plate - 33mm	7 - 12kg
2.7mm TPLO Plate - 43mm	12 - 18kg
3.5mm Mini TPLO Plate - 53mm	18 - 30kg
3.5mm TPLO Plate - 62mm	22 - 40kg
3.5mm Broad TPLO Plate - 78mm	35 - 50kg
3.5mm Jumbo TPLO Plate - 88mm	50kg +